




1 Wash and dry your hands.



2 Pee into the sample tube until it is almost full.



3 Screw the lid on tightly and put the tube in its protective case.




4 Place all your samples and **lab card** into the grey freepost envelope and post in any Royal Mail post box.

What is included in the kit:

This is a simple and reliable test designed to be safely used at home.



Protective case



Sample tube

Contact us:

By text: 07400 370 855

By email: info@fettle.health

When should I take this test?

You should wait 2 weeks after you think you could have been exposed to an infection. This is because Chlamydia and Gonorrhoea can take up to 2 weeks to be detected in a test. You can do this test at any time of day.

Do I need to write my name on the samples?

No, we have already labelled your samples with a unique code, this is to ensure your experience is as easy and discreet as possible.

What happens next?

We will text you when your samples arrive at the lab.

We will then text you your test results as soon as possible (always within 7 days).

If you do have an infection we will advise you how to get treatment and further support.

FAQs:

fettle.health/faqs

POWERED BY
SH:24

Urine sample
Chlamydia + Gonorrhoea

Key information

fettle